

Anti-Pancreasin/Marapsin/PRSS27 antibody



Catalog Number: 100632

Product name

Anti-Pancreasin/Marapsin/PRSS27 antibody

Immunogen

[Human Pancreasin/Marapsin/PRSS27 \(His Tag\) recombinant protein](#)

Specificity

Human Pancreasin / Marapsin / PRSS27

Antibody description

Rabbit monoclonal to
Pancreasin/Marapsin/PRSS27

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Pancreasin / Marapsin / PRSS27 (rh Pancreasin / Marapsin / PRSS27; Q9BQR3; Met1-Lys290).

Formulation

0.2 µm filtered solution in PBS with 5% treha

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1-10µg/ml


IP: 1-4 µg/mg of lysate

Validations



| Lanes | A |
|----------------------------|--|
| Items | |
| Sample (whole cell lysate) | K562 |
| Sample Volume (µg/lane) | 30 |
| Gel | 13% SDS-PAGE reducing gel |
| Recommended Concentration | 1-10µg/ml |
| Secondary Antibody | DyLight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution. |
| Explanation | Developed using Odyssey imaging system. Predicted band size : 32 kDa Observed band size : 32 kDa |

Pancreasin / Marapsin / PRSS27 Antibody, Rabbit



| Lanes | A | B |
|----------------------------|---|------|
| Items | | |
| Sample (whole cell lysate) | Jurkat | Hela |
| Sample quantity | 0.5 mg | |
| IP antibody quantity | 2 µg | |
| Protein G agarose | 15 µl of 50% Protein G Agarose | |
| Gel | 13% SDS-PAGE reducing gel | |
| Primary antibody | ABHD148 antibody at 10 µg/ml [Cat# 13913-R011] | |
| Secondary antibody | DyLight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution. | |

MAb, Western blot

Pancreasin / Marapsin / PRSS27 Antibody, Rabbit
MAb, Immunoprecipitation