

Anti-IL-1R1/CD121a antibody



Catalog Number: 101310

Product name

Anti-IL-1R1/CD121a antibody

Immunogen

[Mouse IL-1R1/CD121a \(His Tag\) recombinant protein](#)

Specificity

Mouse IL-1R1 / CD121a

Antibody description

Rabbit monoclonal to IL-1R1/CD121a

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL-1R1 / CD121a (rM IL-1R1 / CD121a; P13504; Met1-Lys338) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of

water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



IL1R1 Antibody (PE), Rabbit MAb, Flow cytometric analysis

Flow cytometric analysis of Mouse IL1R1(CD121a) expression on NIH/3T3 cells. NIH/3T3 cells were stained with PE-conjugated anti-mouse IL1R1(CD121a). The histogram were derived from gated events with the forward and side light-scatter characteristics of intact cells.