

# Anti-IL10RB / IL-10RB antibody



Catalog Number: 101534

## Product name

Anti-IL10RB / IL-10RB antibody

## Immunogen

[Human IL10RB / IL-10RB \(His Tag\) recombinant protein](#)

## Specificity

Human IL10RB / IL-10RB

## Antibody description

Rabbit polyclonal to IL10RB / IL-10RB

## Preparation

Produced in rabbits immunized with purified, recombinant Human IL10RB / IL-10RB (rh IL10RB / IL-10RB; NP\_000619.3; Met1-Ser220). IL10RB / IL-10RB specific IgG was purified by Human IL10RB / IL-10RB affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB


## Dilutions

WB: 10-20 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human IL10RB / IL-10RB. The detection limit for Human IL10RB / IL-10RB is approximately 0.00245 ng/well.

## Validations



IL10RB / IL-10RB Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot