

Anti-TIMD4 antibody



Catalog Number: 102582

Product name

Anti-TIMD4 antibody

Immunogen

[Mouse TIMD4 \(His Tag\) recombinant protein](#)

Specificity

Mouse TIMD4 / TIM4

Antibody description

Rabbit polyclonal to TIMD4

Preparation

Produced in rabbits immunized with purified, recombinant Mouse TIMD4 (rM TIMD4; NP_848874.3; Met 1-Thr 279). TIMD4 specific IgG was purified by mouse TIMD4 affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 10-30 μ g/mL

ELISA: 0.5-1 μ g/mL

This antibody can be used at 0.5-1 μ g/mL with the appropriate secondary reagents to detect Mouse TIMD4. The detection limit for Mouse TIMD4 is 0.00245 ng/well.

Validations



TIM4 / TIMD4 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot