

Anti-CD226/DNAM-1 antibody



Catalog Number: 103367

Product name

Anti-CD226/DNAM-1 antibody

Immunogen

[Mouse CD226/DNAM-1 \(His Tag\) recombinant protein](#)

Specificity

Mouse CD226 / DNAM-1

Antibody description

Rabbit monoclonal to CD226/DNAM-1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD226 / DNAM-1 (rM CD226 / DNAM-1; NP_848802.2; Met1-Pro254) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be

disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



CD226 / DNAM-1 Antibody (PE), Rabbit MAb, Flow cytometric analysis

Analysis of CD226 expression on spleen lymphocytes.