

# Anti-Arylsulfatase A / ARSA antibody



Catalog Number: 103922

## Product name

Anti-Arylsulfatase A / ARSA antibody

## Immunogen

[Human Arylsulfatase A / ARSA \(His Tag\) recombinant protein](#)

## Specificity

Human Arylsulfatase A / ARSA  
**Has cross-reactivity** in ELISA with Mouse ARSA

## Antibody description

Rabbit monoclonal to Arylsulfatase A / ARSA

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Arylsulfatase A / ARSA (rh Arylsulfatase A / ARSA; NP\_000478.2; Met 1-Ala 507).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 2-10 µg/mL

## Validations



Items	Lanes	A
Sample (whole cell lysate)		HepG2
Sample Volume (µg/lane)		60
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		2-10 µg/ml
Secondary Antibody		Dylight 800 labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.

Arylsulfatase A / ARSA Antibody, Rabbit MAb, Western blot