

Anti-TH / Tyrosine Hydroxylase antibody



Catalog Number: 104090

Product name

Anti-TH / Tyrosine Hydroxylase antibody

Immunogen

[Mouse TH / Tyrosine Hydroxylase \(His Tag\) recombinant protein](#)

Specificity

Mouse TH / Tyrosine Hydroxylase

Antibody description

Rabbit polyclonal to TH / Tyrosine Hydroxylase

Preparation

Produced in rabbits immunized with purified, recombinant Mouse TH / Tyrosine Hydroxylase (rM TH / Tyrosine Hydroxylase; P24529; Pro 2-Ser 498). TH / Tyrosine Hydroxylase specific IgG was purified by Mouse TH / Tyrosine Hydroxylase affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse TH.

IHC-P: 0.1-2 µg/mL

Validations



TH / Tyrosine Hydroxylase Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of mouse TH in mouse brain with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).