

Anti-HSPD1/HSP60 antibody



Catalog Number: 104314

Product name

Anti-HSPD1/HSP60 antibody

Immunogen

[Mouse HSPD1/HSP60 \(His Tag\) recombinant protein](#)

Specificity

Mouse HSPD1 / HSP60

Antibody description

Rabbit polyclonal to HSPD1/HSP60

Preparation

Produced in rabbits immunized with purified, recombinant Mouse HSPD1 / HSP60 (rM HSPD1 / HSP60; NP_034607.3; Leu2-Phe573). Total IgG was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-2 µg/mL

Validations



HSPD1 / HSP60 Antibody, Rabbit PAb, Immunohistochemistry

Immunochemical staining of mouse HSP60 in mouse intestine with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded



sections).

HSPD1 / HSP60 Antibody, Rabbit PAb, Immunohistochemistry

Immunochemical staining of mouse HSP60 in mouse lung with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).