

# Anti-B3GAT3 / GLCATI antibody



Catalog Number: 105134

## Product name

Anti-B3GAT3 / GLCATI antibody

## Immunogen

[Human B3GAT3 / GLCATI \(His Tag\) recombinant protein](#)

## Specificity

Human B3GAT3 / GLCATI

## Antibody description

Rabbit polyclonal to B3GAT3 / GLCATI

## Preparation

Produced in rabbits immunized with purified, recombinant Human B3GAT3 / GLCATI (rh B3GAT3 / GLCATI; NP\_036332.2; Glu72-Val335). B3GAT3 / GLCATI specific IgG was purified by Human B3GAT3 / GLCATI affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC-P


## Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human B3GAT3 / GLCATI.


IHC-P: 0.1-2 µg/mL

## Validations



B3GAT3 / GLCATI Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human B3GAT3 in human brain with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded



sections).

B3GAT3 / GLCATI Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

# Anti-B3GAT3 / GLCAT1 antibody



Catalog Number: 105134

---

Immunochemical staining of human B3GAT3 in human pancreas with rabbit polyclonal antibody

(0.5  $\mu\text{g}/\text{mL}$ , formalin-fixed paraffin embedded sections).