

Anti-STX3 / Syntaxin 3 antibody



Catalog Number: 106029

Product name

Anti-STX3 / Syntaxin 3 antibody

Immunogen

[Human STX3 / Syntaxin 3 \(His Tag\) recombinant protein](#)

Specificity

Human STX3 / Syntaxin 3

Antibody description

Rabbit polyclonal to STX3 / Syntaxin 3

Preparation

Produced in rabbits immunized with purified, recombinant Human STX3 / Syntaxin 3 (rh STX3 / Syntaxin 3; NP_004168.1; Met1-Asn289). STX3 / Syntaxin 3 specific IgG was purified by Human STX3 / Syntaxin 3 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human STX3 / Syntaxin 3.

IHC-P: 0.1-2 µg/mL

Validations



STX3 / Syntaxin 3 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human STX3 in human kidney with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).