

# Anti-STK4/MST1 antibody



Catalog Number: 106139

## Product name

Anti-STK4/MST1 antibody

## Immunogen

[Human STK4/MST1 \(His Tag\) recombinant protein](#)

## Specificity

Human STK4 / MST1

## Antibody description

Rabbit monoclonal to STK4/MST1

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human STK4 / MST1 (rh STK4 / MST1; Q13043-1; Glu2-Phe487).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 10-30 µg/mL

## Validations



STK4 / MST1 Antibody, Rabbit MAb, Immunocytochemistry

Immunocytochemical staining of human STK4 in human stomach carcinoma with rabbit monoclonal antibody (20 µg/mL, formalin-fixed paraffin embedded sections).