

Anti-FLT-3/CD135/FLK-2 antibody



Catalog Number: 106207

Product name

Anti-FLT-3/CD135/FLK-2 antibody

Immunogen

[Human FLT-3/CD135/FLK-2 \(His Tag\) recombinant protein](#)

Specificity

Human FLT3 / CD135 / FLK2

No cross-reactivity in ELISA with E.coli cell lysate

Antibody description

Mouse monoclonal to FLT-3/CD135/FLK-2

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human FLT3 / CD135 / FLK2 (rh FLT3 / CD135 / FLK2; NP_004110.2; Met 1-Asn 541) and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of

receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



CRABP2 Antibody, Mouse MAb, Immunohistochemistry

FLT3 / CD135 / FLK2 Antibody (APC), Mouse MAb, Flow cytometric analysis