

RAB31-specific antibody



Catalog Number: 106966

Product name

RAB31-specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB31-specific Rabbit Polyclonal antibody.
Positive WB detected in human brain tissue.
Positive IHC detected in human heart tissue.
Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006868).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations




human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114430(RAB31-specific antibody) at dilution of



1:500

Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:114430(RAB31-specific Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:114430(RAB31-specific Antibody) at dilution of 1:50 (under 40x lens)