

# HDAC6-specific antibody

Catalog Number: 106979

## Product name

HDAC6-specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HDAC6-specific Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HT-1080 cells. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 130 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001321231). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111381(HDAC6-specific antibody) at dilution of 1:200



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111381(HDAC6-specific Antibody) at dilution of 1:50 (under 40x lens)