

ING4-specific antibody



Catalog Number: 106983

Product name

ING4-specific antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

ING4-specific Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human brain tissue, human liver cancer tissue. Positive IP detected in HeLa cells. Positive WB detected in Jurkat cells, HEK-293 cells, HeLa cells, mouse colon tissue. Observed molecular weight by Western-blot: 29-32 kDa

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111796(ING4-specific antibody) at dilution of 1:500

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001127582). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF, IP


Dilutions

Recommended Dilution:

WB: 1:200-1:2000




IP Result of anti-ING4-specific (IP:Catalog No:111796, 3ug; Detection:Catalog No:111796 1:300) with HeLa cells lysate 2500ug.



ING4-specific antibody

Catalog Number: 106983

Immunohistochemical of paraffin-embedded human brain using Catalog No:111796(ING4-specific antibody) at dilution of 1:50 (under 40x



lens)

Immunofluorescent analysis of HepG2 cells using Catalog No:111796(ING4-specific Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG