

# BAG1L-speicfic antibody



Catalog Number: 106998

---

## Product name

BAG1L-speicfic antibody

## Specificity

Human; other species not tested.

## Antibody description

BAG1L-speicfic Rabbit Polyclonal antibody.  
Positive IF detected in HepG2 cells.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_004323).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of HepG2 cells, using BAG1 antibody Catalog No:108407 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).