

Claudin 9-specific antibody



Catalog Number: 107004

Product name

Claudin 9-specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Claudin 9-specific Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_020982). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109343(CLDN9-specific antibody) at dilution of 1:100



IP Result of anti-CLDN9-specific (IP:Catalog No:109343, 4ug; Detection:Catalog No:109343 1:300) with HeLa cells lysate 880ug.