

BRD8 antibody



Catalog Number: 107022

Product name

BRD8 antibody

Specificity

Human; other species not tested.

Antibody description

BRD8 Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of BRD8 recombinant protein (Accession Number: BC008039). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107022(BRD8 Antibody) at dilution of 1:1000