

BSAP,PAX5 antibody



Catalog Number: 107025

Product name

BSAP,PAX5 antibody

Specificity

Human; other species not tested.

Antibody description

BSAP,PAX5 Mouse Monoclonal antibody. Positive WB detected in Ramos cells, Daudi cells, Raji cells. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of BSAP,PAX5 recombinant protein (Accession Number: NM_016734). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Ramos cells were subjected to SDS PAGE followed by western blot with Catalog No:107025(BSAP,PAX5 Antibody) at dilution of 1:1500