

C1QC antibody

Catalog Number: 107026

Product name

C1QC antibody

Specificity

Human; other species not tested.

Antibody description

C1QC Mouse Monoclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of C1QC recombinant protein (Accession Number: NM_001347619). Purification method: Thiophilic affinity chromatograph.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107026(C1QC Antibody) at dilution of 1:2000