

# C4 (gamma chain) antibody



Catalog Number: 107028

## Product name

C4 (gamma chain) antibody

## Specificity

Human; other species not tested.

## Antibody description

C4 (gamma chain) Mouse Monoclonal antibody.  
Positive WB detected in human plasma tissue.  
Observed molecular weight by Western-blot: 31 kDa

## Preparation

This antibody was obtained by immunization of C4 (gamma chain) recombinant protein (Accession Number: BC012372). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

## Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107028(C4A Antibody) at dilution of 1:8000