

# CD109 antibody



Catalog Number: 107038

## Product name

CD109 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD109 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 162 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM\_005248659). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107038(CD109 Antibody) at dilution of 1:4000