

CD11a/Integrin alpha-L antibody



Catalog Number: 107039

Product name

CD11a/Integrin alpha-L antibody

Specificity

Human; other species not tested.

Antibody description

CD11a/Integrin alpha-L Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of CD11a/Integrin alpha-L recombinant protein (Accession Number: NM_001114380). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (CD11a/Integrin alpha-L Antibody) at dilution of 1:500 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (CD11a/Integrin alpha-L Antibody) at dilution of 1:500 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).