

# CD13 antibody



Catalog Number: 107041

## Product name

CD13 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD13 Mouse Monoclonal antibody. Positive WB detected in U-937 cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 150 kDa

## Preparation

This antibody was obtained by immunization of CD13 recombinant protein (Accession Number: XM\_011521473). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:107041(ANPEP Antibody) at dilution of 1:4000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (ANPEP Antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (ANPEP Antibody) at dilution of 1:100 (under 40x lens)