

CD2 antibody



Catalog Number: 107049

Product name

Monoclonal

CD2 antibody

Ig Type

Specificity

Mouse IgG1

Human; other species not tested.

Applications

Antibody description

ELISA, WB

CD2 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 45~55 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of CD2 recombinant protein (Accession Number: NM_001767). Purification method: Caprylic acid/ammonium sulfate precipitation.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107049(CD2 antibody) at dilution of 1:400