

# Annexin V antibody



Catalog Number: 107058

## Product name

Annexin V antibody

## Specificity

Human; other species not tested.

## Antibody description

Annexin V Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, K-562 cells. Positive FC detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of Annexin V recombinant protein (Accession Number: NM\_001154). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF, FC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107058(ANXA5 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:107058(ANXA5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)



1X10<sup>6</sup> HeLa cells were stained with 0.2ug ANXA5 antibody (Catalog No:107058, red) and

# Annexin V antibody



Catalog Number: 107058

---

control antibody (blue). Fixed with 90% MeOH  
blocked with 3% BSA (30 min). Alexa Fluor 488-

conjugated AffiniPure Goat Anti-Mouse IgG(H+L)  
with dilution 1:1500.