

# APOBEC3G antibody



Catalog Number: 107063

## Product name

APOBEC3G antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

APOBEC3G Mouse Monoclonal antibody. Positive WB detected in A375 cells, human testis tissue, PC-3 cells. Positive IHC detected in human breast cancer tissue. Positive IF detected in C6 cells. Observed molecular weight by Western-blot: 46 kDa

## Preparation

This antibody was obtained by immunization of APOBEC3G recombinant protein (Accession Number: NM\_021822). Purification method: Protein A purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB, IHC, IF

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000


## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:107063(APOBEC3G antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107063(APOBEC3G antibody) at dilution of 1:50 (under 10x lens)




# APOBEC3G antibody



Catalog Number: 107063

---

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107063(APOBEC3G antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of C6 cells using Catalog No:107063(APOBEC3G Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Mouse IgG