

APOBEC3G antibody



Catalog Number: 107063

Product name

APOBEC3G antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

APOBEC3G Mouse Monoclonal antibody. Positive WB detected in A375 cells, human testis tissue, PC-3 cells. Positive IHC detected in human breast cancer tissue. Positive IF detected in C6 cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of APOBEC3G recombinant protein (Accession Number: NM_021822). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:1000


Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:107063(APOBEC3G antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107063(APOBEC3G antibody) at dilution of 1:50 (under 10x lens)




APOBEC3G antibody



Catalog Number: 107063

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107063(APOBEC3G antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of C6 cells using Catalog No:107063(APOBEC3G Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Mouse IgG