

ATG5 antibody



Catalog Number: 107070

Product name

ATG5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ATG5 Mouse Monoclonal antibody. Positive WB detected in K-562 cells, HeLa cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of ATG5 recombinant protein (Accession Number: NM_004849). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107070(ATG5 Antibody) at dilution of 1:1000