

B2M antibody

Catalog Number: 107076

Product name

B2M antibody

Specificity

Human; other species not tested.

Antibody description

B2M Mouse Monoclonal antibody. Positive IHC detected in human prostate cancer tissue, human renal cell carcinoma tissue. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 13 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004048). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107076(B2M Antibody) at dilution of 1:4000



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:107076 (B2M Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:107076 (B2M Antibody) at dilution of 1:200 (under 40x lens)