

# BBS2 antibody

Catalog Number: 107078

## Product name

BBS2 antibody

## Specificity

Human, Dog; other species not tested.

## Antibody description

BBS2 Mouse Monoclonal antibody. Positive WB detected in K-562 cells. Positive IF detected in MDCK cells. Observed molecular weight by Western-blot: 80-90 kDa

## Preparation

This antibody was obtained by immunization of BBS2 recombinant protein (Accession Number: NM\_031885). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG3

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



Immunofluorescent images of MDCK cells stained with BBS2 mouse mAb (Catalog No:107078) and ARL13B rabbit pAb (Catalog No:108196) at dilution of 1:50, further stained with Alexa Fluor 488-conjugated AffiniPure-Goat anti-Mouse IgG(H+L) for Catalog No:107078, and Alexa Fluor 594-conjugated AffiniPure Goat Anti-mouse



IgG(H+L) for Catalog No:108196.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107078(BBS2 Antibody) at dilution of 1:1000