

HRP-conjugated beta Actin antibody



Catalog Number: 107088

Product name

HRP-conjugated beta Actin antibody

Specificity

Human, Mouse, Rat, Zebra Fish; other species not tested.

Antibody description

HRP-conjugated beta Actin Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_001101). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

Validations



Western blot analysis of beta actin in various tissues and cell lines using Proteintech antibody HRP-60008 at a dilution of 1:5000.