

# AEBP2 antibody



Catalog Number: 107101

## Product name

AEBP2 antibody

## Immunogen

[Human AEBP2 Recombinant protein \(His tag\)](#)

## Specificity

human; other species not tested.

## Antibody description

AEBP2 Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of AEBP2 recombinant protein (Accession Number: BC015624). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107561(AEBP2 antibody) at dilution of 1:1000