

AIRE antibody



Catalog Number: 107104

Product name

AIRE antibody

Specificity

Human, Pig; other species not tested.

Antibody description

AIRE Mouse Monoclonal antibody. Positive WB detected in human spleen tissue, pig spleen tissue. Observed molecular weight by Western-blot: 56-59 kDa

Preparation

This antibody was obtained by immunization of AIRE recombinant protein (Accession Number: NM_000658). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:107564(AIRE Antibody) at dilution of 1:2000