

Product name

Monoclonal

AIRE antibody

Ig Type**Specificity**

Mouse IgG1

Human, Pig; other species not tested.

Applications**Antibody description**

ELISA, WB

AIRE Mouse Monoclonal antibody. Positive WB detected in human spleen tissue, pig spleen tissue. Observed molecular weight by Western-blot: 56-59 kDa


Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of AIRE recombinant protein (Accession Number: NM_000658). Purification method: Protein G purified.

Validations**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:107564(AIRE Antibody) at dilution of 1:2000

Clonality