

# CD3 gamma antibody



Catalog Number: 107130

## Product name

CD3 gamma antibody

## Specificity

Human; other species not tested.

## Antibody description

CD3 gamma Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 20-25 kDa

## Preparation

This antibody was obtained by immunization of CD3 gamma recombinant protein (Accession Number: NM\_000073). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:50-1:500


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107130(CD3G Antibody) at dilution of 1:4000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107130(CD3G Antibody) at dilution of 1:400 (under 10x lens). heat mediated antigen retrieved with Sodium Citrate buffer (pH6).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107130(CD3G Antibody) at dilution of 1:400 (under 40x lens). heat mediated antigen retrieved

# CD3 gamma antibody



Catalog Number: 107130

---

with Sodium Citrate buffer (pH6).