

CD47 antibody



Catalog Number: 107138

Product name

CD47 antibody

Specificity

Human; other species not tested.

Antibody description

CD47 Mouse Monoclonal antibody. Positive WB detected in SKOV-3 cells, HT-1080 cells, HUVEC cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of CD47 recombinant protein (Accession Number: NM_198793). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107138(CD47 Antibody) at dilution of 1:3000