

# CD7 antibody

Catalog Number: 107141

---

## Product name

CD7 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CD7 Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue.

## Preparation

This antibody was obtained by immunization of CD7 recombinant protein (Accession Number: NM\_006137). Purification method: .

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:100-1:400

## Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107141 (CD7 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer (pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107141 (CD7 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer (pH9).