

CDK6 antibody



Catalog Number: 107148

Product name

CDK6 antibody

Specificity

Human; other species not tested.

Antibody description

CDK6 Mouse Monoclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in Jurkat cells, THP-1 cells. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of CDK6 recombinant protein (Accession Number: BC027989). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107148(CDK6 Antibody) at dilution of 1:2000



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Catalog No:107148(CDK6 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)