

CDX2 antibody



Catalog Number: 107149

Product name

Monoclonal

CDX2 antibody

Ig Type

Specificity

Mouse IgG1

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

CDX2 Mouse Monoclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 33-40KD

Dilutions


Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of CDX2 recombinant protein (Accession Number: NM_001265). Purification method: Protein G purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107149(CDX2 antibody) at dilution of 1:1000

Clonality