

CHRNA7 antibody



Catalog Number: 107156

Product name

CHRNA7 antibody

Immunogen

[Human CHRNA7 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

CHRNA7 Mouse Monoclonal antibody. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of CHRNA7 recombinant protein (Accession Number: XM_017021884). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107156(CHRNA7 antibody) at dilution of 1:1000