

# CLEC4G antibody



Catalog Number: 107161

## Product name

CLEC4G antibody

## Specificity

Human; other species not tested.

## Antibody description

CLEC4G Mouse Monoclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of CLEC4G recombinant protein (Accession Number: NM\_198492). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107161(CLEC4G Antibody) at dilution of 1:1000