

Complement factor B antibody



Catalog Number: 107169

Product name

Complement factor B antibody

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

Complement factor B Mouse Monoclonal antibody.
Positive IHC detected in human liver tissue.
Positive WB detected in human blood tissue.
Observed molecular weight by Western-blot: 93 kDa




Preparation

This antibody was obtained by immunization of Complement factor B recombinant protein (Accession Number: NM_001710). Purification method: Thiophilic Affinity Chromatography.

human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:107169(CFB antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Monoclonal

Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:107169(CFB Antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Mouse IgM



Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:107169(CFB Antibody) at dilution of 1:50

Complement factor B antibody



Catalog Number: 107169

(under 40x lens)