

CPT1C antibody



Catalog Number: 107174

Product name

CPT1C antibody

Specificity

human; other species not tested.

Antibody description

CPT1C Mouse Monoclonal antibody. Positive IHC detected in human testis. Positive WB detected in Human brain. Observed molecular weight by Western-blot: 82kD

Preparation

This antibody was obtained by immunization of CPT1C recombinant protein (Accession Number: XM_047438159). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107174(CPT1C antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human testis using Catalog No:107174(CPT1C antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:107174(CPT1C antibody) at dilution of 1:50 (under 40x lens)