

CXCL8/IL8 antibody



Catalog Number: 107185

Product name

CXCL8/IL8 antibody

Specificity

Human; other species not tested.

Antibody description

CXCL8/IL8 Mouse Monoclonal antibody. Positive WB detected in ag10552. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of CXCL8/IL8 recombinant protein (Accession Number: NM_000584). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107185 at dilution of 1:1000