

FABP5 antibody



Catalog Number: 107198

Product name

FABP5 antibody

Specificity

Human; other species not tested.

Antibody description

FABP5 Mouse Monoclonal antibody. Positive WB detected in A375 cells, fetal human brain tissue. Observed molecular weight by Western-blot: 15 kDa

Preparation

This antibody was obtained by immunization of FABP5 recombinant protein (Accession Number: NM_001444). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:107198(FABP5 Antibody) at dilution of 1:8000