

Factor XII antibody



Catalog Number: 107199

Product name

Factor XII antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

[Human F12 Recombinant protein \(His tag\)](#)

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

Factor XII Mouse Monoclonal antibody. Positive IHC detected in human liver tissue. Positive IF detected in HepG2 cells. Positive WB detected in K562 cells. Observed molecular weight by Western-blot: 28 kDa




Preparation

This antibody was obtained by immunization of Factor XII recombinant protein (Accession Number: BC012390). Purification method: Protein A purified.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107199(F12 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human liver using Catalog No:107199(F12 antibody) at dilution of 1:50 (under 10x lens)

Clonality


Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IF, IHC




Dilutions

Factor XII antibody

Catalog Number: 107199

Immunohistochemical of paraffin-embedded human liver using Catalog No:107199(F12 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using F12 antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).