

Factor XIIIa antibody



Catalog Number: 107200

Product name

Factor XIIIa antibody

Specificity

Human, Pig, Mouse; other species not tested.

Antibody description

Factor XIIIa Mouse Monoclonal antibody. Positive WB detected in human placenta tissue, mouse liver tissue, pig liver tissue. Observed molecular weight by Western-blot: 75-80 kDa

Preparation

This antibody was obtained by immunization of Factor XIIIa recombinant protein (Accession Number: NM_000129). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:107200 (Factor XIIIa Antibody) at dilution of 1:4000