

Factor XIIIa antibody



Catalog Number: 107200

Product name

Factor XIIIa antibody

Monoclonal

Ig Type

Mouse IgG2b

Specificity

Human, Pig, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

Factor XIIIa Mouse Monoclonal antibody. Positive WB detected in human placenta tissue, mouse liver tissue, pig liver tissue. Observed molecular weight by Western-blot: 75-80 kDa

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of Factor XIIIa recombinant protein (Accession Number: NM_000129). Purification method: Protein A purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:107200 (Factor XIIIa Antibody) at dilution of 1:4000

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality