

# FAM3C antibody



Catalog Number: 107203

## Product name

FAM3C antibody

## Specificity

Human; other species not tested.

## Antibody description

FAM3C Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 25 kDa

## Preparation

This antibody was obtained by immunization of FAM3C recombinant protein (Accession Number: XM\_047419774). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107203(FAM3C Antibody) at dilution of 1:1000